

PURCHASE ORDER FORM

URODYNAMICS CONSUMABLES and SUPPLIES

ABN 35082177231

Kirkwall Australia Pty Ltd trading as:



Tel: 02 9450 0855

Fax: 02 9450 0844

P.O. Box 290 PYMBLE NSW 2073

7/5 Narabang Way, Belrose NSW 2085

Purchase Order No.		ABN	
Invoice To:		Ship To:	
Suburb		Suburb	
State	P/Code	State	P/Code
Tel:	Authorised by:		
Fax:	Signed (see note 3)		

Code	Units/Box	Description	Qty Req'd	Price Excl GST	Total
Catheters					
601.392	10/Box	5200 Filling Catheter 10F 40cm. Connect to std giving set			
601.393	10/Box	5212 Filling Catheter 12F 40cm. Connect to std giving set			
601.991	25/Box	5145 Abdominal Pressure Catheter 4.5F (direct)			
601.389	25/Box	5400 Abdominal Pressure Catheter 4.5F (Sheath)			
601.391	25/Box	5100 Bladder Pressure catheter 4.5Fr			
601.657	10/Box	5706 Dual-bladder/urethral-Press catheter 8F			
601.515	10/Box	5716 Dual-bladder Fill & Press catheter 6F			
601.427	10/Box	5331 Triple Lumen Filling Plus Bladder & Urethral pressures			
601.731	25/Box	5415 Rectal Balloon Catheter 4.5 Fr - non sterile			
Perfusion Sets					
601.516	5/Box	PL-2 Perfusion Set - 2 lines with 2 taps- Neopress			
601.506	5/Box	PL-3 Perfusion Set - 3 lines with 3 taps - Neopress			
Transducers & Domes					
601.676	50/Box	MX960XX (M-M) Pressure transducer Domes			
600.044	50/Box	MX848 (M-F) Pressure transducer Domes For MX860 transducer (Needs MX493 Adapter)			
600.045	50/Box	MX858 (M-M) Pressure transducer Domes (for MX860)			
600.482	50/Box	MX534-1L Stopcock - four way (C kit)			
601.424	50/Box	MX 493 Adapter (M-M) - Luer, sterile (C kit), If using MX848 Domes.			
For Australian Orders: Freight & GST will be added				TOTAL A\$	

NEOPRESS SYSTEM RECOMMENDATION

SINGLE USE CATHETERS

Note 1: SESSION ITEMS typically do not require replacement **between patient tests within a testing session**. These items include Transducer Domes and Perfusion Sets with integral taps. Use of the Neopress Perfusion System simplifies diagnostic equipment set-up and patient testing procedures thus providing for an extremely cost effective patient study solution.

Note 2: Your signature on this document implies that you have on behalf of the END USER read and accepted the Neomedix Consumables Terms and Conditions of Sale and Notice of Disclaimer available at www.neomedix.com.au or available upon request.